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

NO.	Category	Quantity	Note
1	MCI-1 Oxygen Analyzer	One Set	
2	Oxygen Connecting Tube	One Piece	
3	User Manual	One Piece	

● Oxygen Purity Calibration Procedure:

Turn on. Switch to the purity/flow measurement mode. Press and hold the "←" and "M" buttons until the calibration mode selection interface is indicated. (see chart)


Then perform oxygen purity calibration. Calibration modes are divided to two modes, low oxygen calibration and high oxygen calibration. Under calibration, only one mode can be selected to perform calibration.

Low oxygen calibration: use air source for low oxygen point calibration. Ensure that there is no oxygen residue in the MCI-1, when calibration mode selection interface is indicated, press the "M" button to perform the low oxygen calibration, and the screen will returns to purity/flow measurement interface itself.

High oxygen calibration: use oxygen source for high oxygen point calibration. When calibration mode selection interface is indicated, press the "←" button, the MCI-1 will return to the purity/flow display interface. Connect the oxygen (known concentration) tubing to the purity/flow sample inlet of the MCI-1. Then press the "←" button or "M" button and adjust the indicated purity to keep it consistent with the one of the calibrated gas. High oxygen calibration finishes. Then measure other gas.

User Manual For MCI-1 Oxygen Analyzer

Technical Date

Measurement Parameters:

Measurement Item	Measurement Range	Measurement Resolution	Measurement Accuracy	Remark
Oxygen Purity	30%-96%	0.1%	±3%	For gas at the outlet of PSA oxygen concentrator only
Flow	1-10L/min	0.1L/min	±0.2L/min (1-2L)	The indicated measured value is surface pressure
			±0.3L/min (2-10L)	
Pressure	0-3bar	1 kPa	±2%FS	is surface pressure
		0.1bar		
		0.1psi		
Power Supply		One 9V battery		
Measurement Current		33mA		
Standby Current		13µA		
Operating Temperature		5~35°C		
Storage Temperature		-15~60°C		
Dimensions		145x66x48(mm)		
Weight		115(g)		

Malfunction and Troubleshooting

Malfunction	Troubleshooting Method
The display value has obvious error, and the battery icon appears at the top right of the screen	Battery low and replace battery
The flow measurement value is significantly lower than the actual value	Check whether there is air leakage in pipes or tubing
The pressure measured value is significantly lower than the actual value	Check whether there is air leakage in pipes or tubing
No response after connecting measured gas	Check whether the purity/flow sample inlet or pressure sample inlet is connected incorrectly
The purity measurement error is large	Re-calibrate
The indoor measurement error is large in winter	The If the temperature of the environment changes too much, such as entering the room from the outdoors in winter, please put oxygen analyzer in room for more than 15 minutes before use
Unable to boot	Check whether the anode and cathode of the 9V battery are connected well

Warranty

The warranty period of the whole machine: 1 year
During warranty period, we will only provide paid maintenance instead of free maintenance for the following malfunction caused by users.

- Malfunction caused by unauthorized dismantling or refitting;
- Malfunction caused by the careless drop in use of handling;
- Malfunction caused by lack of reasonable maintenance;
- Malfunction caused by failure to operate the product in accordance with the proper guide in this user manual.

Product Packing List (including list of replaceable parts)

NO.	Category	Quantity	Note
1	MCI-1 Oxygen Analyzer	One Set	
2	Oxygen Connecting Tube	One Piece	
3	9V Battery	One Piece	Replaceable parts
4	User Manual	One Piece	
5	Certificate of Qualification	One Piece	

SHENYANG AERTI TECH CO.,LTD.
No.77-1, 13th Road, Shenyang Economic & Technological Development Area, Shenyang City, 110027
Liaoning, P.R.China
Tel: +86-024-31229971-8002
Email: anna@aerti.com.cn

IMPORTED BY

MOBIAC S.A
KAOIANA AKPOTHPIOY, XANIA - KPHTH
KATHIANA AKROTIRIOU CHANIA - CRETE
T: +30 28210 63222 | WEB: www.mobiaccare.com

Contents

Product introduction.....	2
● Intended Use.....	3
● Component Identification and Description	3
Notes.....	4
Operating Instructions.....	4
● ON/OFF.....	4
● Mode Switch.....	5
● Oxygen Purity/Flow Measurement.....	5
● Pressure Measurement.....	5
● Oxygen Purity Calibration Procedure.....	6
Technical Date.....	7
Malfunction and Troubleshooting.....	7

Product Introduction

MCI-1 oxygen analyzer is an measurement instrument developed by SHENYANG AERTI TECH CO.,LTD for PSA oxygen concentrator. It is designed to measure oxygen concentration, flow and outlet pressure of oxygen concentrators. The MCI-1 provides unparalleled performance and reliability from its advanced design that includes the following features and operational benefits:


- *Long battery with 9V *Portable and smart
- *Convenient to measure *Accurate oxygen measurement

2

Intended Use

MCI-1 is for use only by trained personnel.
MCI-1 Oxygen Analyzer is a tool only used to measure output gas at the outlet of a medical molecular sieve oxygen concentrator (oxygen purity, flow and pressure)
MCI-1 Oxygen Analyzer is not intended to be used by patients who are prescribed oxygen, nor is it intended to continuously monitor or confirm oxygen delivery to a patient.
MCI-1 Oxygen Analyzer is intended to be used in an environment where oxygen concentrators are being serviced or repaired.

Component Identification and and Description

3

Notes

- Do not use this product without first completely reading and understanding this user manual.
- Under high oxygen calibration mode, do not use any source other than medical molecular sieve oxygen concentrator. Concentration 85 ~ 93% and flow 2-3 L/min are recommended.
- Do not use this product to measure any gas other than the one at the outlet of medical molecular sieve oxygen concentrator. When measuring purity and flow, gas flow is less than 10L/min, gas pressure is lower than 0.5bar, gas temperature range is 5~35°C. Gas must not be humidified before measured.
- Do not use this product in an open fire or in a flammable environment.
- Not suitable for use in the presence of flammable anesthesia mixture.
- Do not smoke while using the MCI-1 for checking oxygen concentrators.
- Do not use in an MRI environment.
- When not in use for periods greater than 30 days remove the battery to protect the MCI-1 from potential battery leakage.
- Clean the MCI-1 with the soft moist cloth and mild hand or dish soap. Do not drop or immerse the MCI-1 in the water.

4

Operating Instructions

- ON/OFF:
Connected with the battery, default state of the MCI-1 is on, then press the key "⏻" and power off. Or wait about 3 minutes and it will switch itself off. When the MCI-1 is off, press the key "⏻" and it switches on.
- Mode Switch:
When it switches on, the default display is purity/ flow measurement interface. Press the "M" button to switch to the pressure measurement interface. Measurement modes are switched between purity/low measurement mode and pressure measurement mode by pressing the "M" button.
- Oxygen Purity/Flow Measurement:
Turn on. Switch to the purity/flow measurement mode. Connect the gas sample tubing to the purity/flow sample inlet. Allow the oxygen reading to stabilize for approximately 10 seconds before reading the oxygen concentration and flow (see chart)
- Pressure Measurement:
Turn on. Switch to the pressure measurement mode. Connect the gas sample tubing to the pressure sample inlet. Allow the oxygen reading to stabilize for approximately 10 seconds before reading the gas pressure (see chart)

O2 Purity (%)	93.0
Flow (L/min)	3.0

Pressure	
0	kPa
0.0	bar
0.0	psi

5