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CE

## Preface

Please read the user manual carefully before taking the product into use.

- ✧ This manual contains operations, assembly methods, and simple faults solutions.
- ✧ This manual applies to our models: 0811315 / 0811317
- ✧ This manual contains wheelchair maintenance and self-checking methods, please put it in proper place.
- ✧ Please provide this manual for reference when other people are going to use this wheelchair.
- ✧ The annotations and illustrations in this manual might be slightly different with the real parts due to quality improvement or changing design. Please in kind prevail.
- ✧ Contact with your dealer if there is any ambiguity or question.
- ✧ Improper use of any vehicle may lead to injury. Unsafe driving could harm yourselves and others.
- ✧ The electric wheelchair is intended to comfortably transport persons with walking difficulties or no walking abilities.
- ✧ This electric wheelchair is designed to transport 1 person only.

Symbols used in this manual.

Warning symbol

Follow the instructions next to this symbol closely.

Not paying careful attention to these instructions could result in physical injury or damage to the wheelchair or the environment.


# 1. Product Introduction

## 1.2 Product Images

REF:0811317

REF:0811315

### 1.3 Product Information

0811317	Requirements
Frame	Aluminum- Folding (once battery is removed)
Controller	Right or Left Adjustable Controller with lighting system
Motor	320W x 24 V x 2
Battery	2 pieces x 12V 40Ah or 35Ah
Charger	24V 5A
Max Speed	8-10 km/h
Autonomy	30-36 km
Gradeability	10°
Front Wheel	9" x 2" (230x50mm) PU
Rear Wheel	13,5" PU
Armrest	Detachable- Height Adjustable
Footrest	Detachable, with Adjustable Length with safety belt
Seat	Removable Cushion
Brake	Electromagnetic Brakes
Other features	Mudguard, seat belt, anti-tilt casters, Adjustable Head support, Adjustable Backrest

46cm	44cm	50cm	40cm	65cm	110cm	100cm	55cm	13,5" 9'	130kg	67kg

0811315	Requirements
Frame	Steel- Folding
Controller	Right or Left Adjustable Controller
Motor	250W x 24 V x 2
Battery	2 pieces x 12V 20AH
Charger	24V 3A
Max Speed	7 Km/h
Autonomy	15-20 km
Gradeability	10°
Front Wheel	8" x 2" PU (200x50mm)
Rear Wheel	12,5" PU
Armrest	Folding with PU cushion
Footrest	Detachable with Adjustable Length with safety belt
Seat	Removable Cushion
Brake	Electromagnetic Brakes
Other features	Seat belt, anti-tilt casters, Tilting Backrest

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### **Product Performance**

This series of electric wheelchairs are powered by battery, driven by motor. Users control direction and adjust speed by joystick controller. It is suitable for application of low speed, good road condition and small slope.

### **Main Structure**

This series of electric wheelchairs are mainly composed of two parts:

- 1) Frame Part: frame, wheel, footrest, armrest, backrest etc.
- 2) Electrical Part: controller, motor, battery.

### **Application Range**

Our electric wheelchair is for disabled and elderly people using as walking vehicle, is suitable for outdoor middle distance use. Do not run on the motorways.

## **2. Safety Instructions**

### **2.1 Driving Notice**

#### **General Driving Cautions**

- ♦ Please keep your hands on the armrest to control the wheelchair.
- ♦ Please practice in parks or other safe open places until you can use the wheelchair skillfully.
- ♦ Fully practice driving in a safe place, to master principle of moving forward, stopping and turning circles.
- ♦ Before driving to the road, please be accompanied by caregivers and confirm it is safe.
- ♦ Please strictly follow traffic rules of the pedestrian, do not regard yourself as vehicle driver.
- ♦ Please drive on the sidewalk and the zebra crossing. Do not drive on vehicle lane.
- ♦ Steadily drive, to avoid driving in 'Z' line or sharp turn.
- ♦ Please keep pressure of pneumatic tire properly. The abnormal tire pressure may cause unsteady driving or excessive lost of current.

#### **Accompanied by Caregivers or Avoid Driving in These Conditions**

- ♦ Drive in bad weather, such as rainy day, heavy fog, strong wind, snow.
- ♦ In case of wheelchair is wet out, wipe the water immediately.
- ♦ Drive on the bad road condition, such as muddy, trail, sand road, gravel, etc.
- ♦ Drive on crowded road.
- ♦ Drive on no fence side ditch, pond, etc.
- ♦ Across the railway.
- ♦ When you have to across railway, pause at "turn off" to confirm it is safe, and make sure that the tires will not be stuck by railway.
- ♦ The electric wheelchair is only for personal use, do not carry people or goods, and do not for a traction purposes.

#### **Precautions for Uphill and Downhill**

- ♦ Avoid driving on steep slopes. The tilt range should be up to 10 degrees.
- ♦ Keep moving forward when going uphill or downhill.
- ♦ Keep speed low when going downhill.
- ♦ Driving on steps is prohibited.
- ♦ Avoid potholes.

#### **Precautions for the Caregivers**

- ♦ Caregivers should confirm that the user's feet are on the footrests proper position and ensure clothes do not attach to the wheels.
- ♦ Caregivers should push the wheelchair to move forward to keep safe on steep slope or long slope.

## 2.2 Other Notice

### Repair and Reform

If the repair and modification is necessary, please contact your dealer. Do not modify by yourself, it may lead to accident or malfunction.

### Protect Environment

In order to protect environment, do not casually discard the abandoned wheelchair and scrapped battery. Please contact your local dealer for replacement of new battery.

### Keep the Storage Environment Dry

Do not storage wheelchair in damp condition, such as near-by bathroom. If wheelchair got wet in the rain, wipe the water immediately.

*Do not replace parts or material randomly.*

*Do not add weight randomly, in case to cause imbalance.*

*Do not use other vehicle to pull or push the wheelchair when somebody sits on the wheelchair or the wheelchair is in manual mode. When there is abnormal sound, contact service immediately.*

### Protection from Electromagnetic Interference

The wheelchair should away from strong magnetic fields and large inductive electrical equipment. Pay attention to that if there are sources of electromagnetic interference nearby, as far as away from those sources to avoid electromagnetic interference

## 3 Adjustment of Wheelchair

### About item 0811317:

Open the box.

Remove the back cushion, seat cushion, headrest, accessory box.

Turn the handle marked in red in the picture to the left so that the wheelchair can be moved.

Take the electric wheelchair out of the box.

Open the wheelchair.

Adjust the seat back.

Adjust the back straps as shown in the images.

Install the back cushion and seat cushion.

Install the controller in the preferable position.

Open the parts box and place the headrest tube on the frame.

Adjust the position of the headrest.

Install the footrests.

**About item 0811317:**

First step: Take off the battery and unplug the red plug (picture 1-a);

Second step: Lift up the cushion strongly by hands, fold wheelchair (picture 1-b).

a

Picture 1

b

**Connect Wire Line**

Connect the battery box's plug to controller's plug (picture 2).

2

### Adjustment of the controller

Adjust the controller to your ideal position by loosening the nut counterclockwise. Tighten the nut normally when you have completed the adjustment (Picture 3).

Picture 3

### Assemble and Disassemble the Footrest

Install the footrest into the frame tube, and move from outside to medial to the fixed position; For disassemble, push the plastic part to outside then lift up footrest (picture 4).

Picture 4

### Adjustment of Castor

Adjust castor angle by unscrew the nut (picture 5).

Picture 5

### Switch between Manual Mode and Electric Mode

 *When changing from manual mode to electric mode, please confirm both side clutches are inside the grooves. Any side clutch is not fully into groove, may cause danger of rollover.*

Manual mode: turn off power and stop the wheelchair, spin the adjusting levers, make it point to the yellow label. (Picture 6-1)

Electric mode: spin the adjusting levers, and make it point to the green label. (Picture 6-2)


Picture 6-1

Picture 6-2

### **Adjusting the backrest angle**

Pulling the strap "1" releases the pins "2" on the right and left, and adjusts the inclination of the back to the positions "3" as shown in the image below.

The inclination of the seat back is adjustable from 100° to 130°, with 4 adjustment positions. The seat back tilt is 100° when the pin is in the position shown in the image below.

### **About item 0811315:**

For wheelchair 0811315 you do not have to remove the batteries to collect it, take out the seat cushion and pull the electric wheelchair to close it

## **4 Battery and Charger**

### **4.1 How to fix the battery box.**

#### **About item 0811315:**

The batteries are mounted on the wheelchair. It can be folded and stored without removing the batteries


### About item 0811317

STEP ONE: put the battery in the battery box and fix the bolts

For "A" fix the bolts like this 3. For "B"fix the bolts like this (upside down)

STEP TWO: fix the battery box into the frame.

## 4.2 Charger Requirements

The charger is used for battery charging. During electric wheelchair is charging, do not use it.

### Charger Technical Data

**0811315:** Input Voltage: AC220V  $\pm$ 10% Output Voltage: 24V/DC Output Current: 3A

**0811317:** Input Voltage: AC220V  $\pm$ 10% Output: 24V/DC Output Current: 5A

## 4.3 Usage of Charger

In order to charge the battery, connect the charger with power supply and battery box's plug. Do as following instructions to complete charging process:

*Step 1:* Make sure the charger groove is not blocked.

*Step 2:* Make sure the electric wheelchair is power off.

*Step 3:* Unplug the plugs that connect the battery box and controller.

*Step 4:* Connect the charger's output plug to the battery box's power plug.

*Step 5:* Connect the charger's main plug to the power supply and the red light will be alight. Fully charging needs 8 to 10 hours, do not overcharge more than 24 hours.

*Step 6:* As to prevent cut short battery life, please charge the battery at least once per month when the wheelchair is not being used.

*Do not stop charging until charge process is finished. Repeat using the battery, which is not fully charged, will shorten the battery life, so the battery should be fully charged as much as possible. When battery is fully charged, the power indicator will turn to green. Do not stop charging before fully charged.*

*After finish charging, turn off the power supply, otherwise the battery will discharge slowly. Do not charge more than 24 hours. Over-charging is dangerous.*

*Users should follow the following rules to avoid charging dangers:*

*The electric wheelchair does not include the charger, please use national standard charger that the output voltage is 24V, 3A.*

*It should be well ventilated when charging. Do not expose the wheelchair to sunlight and humid environment.*

*The charging environment temperature range is from 10°C to 50°C. If it is out of the environment temperature range, the battery is unable to function well, and can easily make battery damaged.*

*It is normal for the fan making a sound during charging. It is for cooling the charger, please do not worry about it.*

*Prevent liquid go into the charger during charging. Moreover, do not place charger on the flammable items, such as fuel, footrest or seat cushion.*

#### **4.4 The Usage and Maintenance of Battery**

Wrong operation of replacing battery may cause danger of explosion. Only the same or recommended type of battery is suitable for replacement. In addition, please make sure the battery poles are correct. Key points for prolong the battery life: charge frequently, to keep battery power full. It is better to fully charge the battery if the wheelchair is not being used. If stop using for a long time, it is better to charge twice per month.

### **5. System Diagnose**

When the indicator LED lights are blinking, it means the wheelchair has abnormalities. The abnormalities may occur in following parts: motors, brakes, battery, wire connections, etc.

Through the product's inner information consultation, the property of the abnormal condition can be detected by the diagnosis signal. The abnormal situation can be detected without other service tools.

#### **Audio Signal Indication**


Description of LED light	The meaning of LED light	Explanation & Solution
All LED lights are unlit with no sound	The power is off, wheelchair in standby or sleep mode. Power is poor contacted. Fuse is tripped or burned out.	
All LED lights are lit	The power is turned on, and self diagnose is passed, the electric wheelchair can work well.	Less LED lights lit, less battery power remaining.
The leftmost red LED light is lit	The battery power is extremely inadequate.	To charge immediately. Or battery is malfunction, and is not able to charge.
Two short beeps with five indicator lights blink twice	The left side motor is malfunction.	The left side motor is poor connected or the wire is disconnected.
Four short beeps with five indicator lights blink four times	The right side motor is malfunction.	The right side motor is poor connected or the wire is disconnected.

Six short beeps with five indicator lights blink six times	Controller is in over-current protection status.	Check the brakes, and check if the motor drive mechanism is stuck. Check the current by ammeter, if it is not excessive current, maybe the controller is malfunction.
Seven short beeps with five indicator lights blink seven times	Joystick is malfunction	Joystick does not reset, or the connector is loose.
Eight short beeps with five indicator lights blink eight times	Controller is malfunction.	Please consult your dealer for maintenance.
Nine short beeps with five indicator lights blink nine times	Controller is malfunction.	Please consult service center for maintenance.

## 6. Controller

### 6.1 Controller Panel

### 6.2 Controller Usage

#### Power Switch

Press this button, the battery power gauge lights will turn on from left to right. Press again, all the LED lights are off.

*In some emergency, you can directly turn off power by press power button.*

### **Sleep Mode**

If the joystick has no operation more than 20 minutes, the power turn off automatically, and the system in the sleep mode, system will be woken up from sleep mode by press power button.

### **Speed adjustment**

According to user's habits and the circumstances, the wheelchair driving speed is adjustable.

Adjust speed by press decrease button or increase button.

Speed is divided into five sections, ranged by 20%, 40%, 60%, 80%, and 100% of the max speed.

First gear of speed: 20% of the max speed.

Fifth gear of speed: 100% of the max speed.

### **Horn Button**

Press the horn button, the horn will sound until you release the button.

### **Battery Power Gauge**

After turn on the power, the battery power gauge is working. The battery power gauge also shows the battery remaining power capacity.

As the picture shown, the battery power is full.

When there is only the red or yellow LED light lit, the battery should be charged. And for a long distance driving, the battery should be fully charged. If only the red LED light is lit, the battery power is extremely inadequate, users need to charge battery as soon as possible.

### **Use of joystick**

The wheelchair moving direction is controlled by joystick. The moving extent of joystick also controls the moving speed.

*Please turn on or turn off the power when the joystick is in the middle position. Otherwise, the controller will let out wrong operation signal.*

*Release the joystick to middle and this mistake will stop. If the mistake is continuing, perhaps this part is malfunction, do not use it and contact with your dealer for maintenance.*

### **Directional light buttons**

The buttons are for selecting right/left “flash”

## **7. Maintenance**

Maintenance of wheelchair includes clean wheelchair, check wheel and battery, and charge battery. For further maintenance, please contact with your dealer. Our suggestion is checking the wheelchair biannual every year.

### **7.1 Clean and Battery Maintenance**

#### **Clean Wheelchair**

Please clean the wheelchair regularly.

Clean the parts, which are frequently touched with user's body (such as set cushion, armrest, controller) by a clean and slight wet cloth. Do not use organic solvents to clean.

If it is a patient user, the wheelchair should be cleaned once per week. If the wheelchair is used by an infectious victim, it should be cleaned and disinfected by disinfectants.

#### **Wheels**

Check the tires' air pressure and their wear condition regularly. When the tire tread pattern depth reduces to 1mm, please replace new tire.

#### **Battery**

To make sure that the battery is fully charged frequently. In order to prolong battery life, we strongly recommend users do not to charge until battery run out of power total.

### **7.2 Wheelchair Malfunction and Checking**

When malfunctions happen during the wheelchair is working, please turn off the power before check.

Symptom: completely lose power, and all the LED lights on the controller panel are off.

#### **Check Step:**

*Step 1:* To check if the controller's plug is loose.

*Step 2:* To check if the connection of controller's plug and battery box is loose. Please reinsert the plug connector (hold the plug when pull out the plug. Do not pull the wire to avoid unnecessary damage to the wire line). After above checking, if the wheelchair is still not able to restore the power, or if users have any question for above checking, please connect with your dealer.

### 7.3 Maintenance Checking

The following is a checking list. Electric wheelchair should be ordinal checked according to our suggestion. When you get on or get off the wheelchair, some self checking is done automatically. For more your attention, we particularly list these self-checking items in A area.

<b>A area</b>	Before use, please check if following parts are correct: 1.Backrest 2. Armrest 3. Controller position 4.Footrest 5. battery power 6.Clutch/adjusting lever for conversion between manual mode and electric mode
<b>B area</b>	Check the following parts monthly, to avoid original parts loosing or wearing 1. Screws 2. Brakes 3. Clutch/adjusting lever for conversion between manual mode and electric mode 4. The front and rear wheels and their tread pattern depth 5. Connectors' of controller and charger
<b>C area</b>	For safety, semiannual overall maintenance is needed.

### 8. Condition for Transport and Storage

During transport and storage, the electric wheelchair should be correctly placed as labeled indication.

- ♦ Transportation should be avoided moisture and sunlight and away from heat resource.
- ♦ Storage condition:
  - Environment temperature  $-40^{\circ}\text{C} \sim +55^{\circ}\text{C}$ ;
  - Relative humidity  $\leq 80\%$ ;
  - Air Pressure  $86\text{kPa} \sim 106\text{kPa}$ .

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## 10. Warranty

The product is guaranteed for a period of two\* (2) years from the purchasing date.

The warranty covers factory defects & **DOES NOT** cover damages resulting from misuse, poor maintenance, modification, overuse or non-compliance with the product's operating instructions. Also **DOES NOT** cover the parts that can wear out during use or over time like brakes, wheels and plastic parts.

Damages or defects caused by: natural disasters, unauthorized maintenance or repair, supply problems, carriers are not covered by the warranty.

**\*SIX (6) months for the battery**

No service or spare part is covered unless it has been inspected by an authorized service department or the manufacturer of the product. Shipping costs for goods of warranty are to be paid by the user- customer or by dealer

### Warranty Claim Form

<b>BUYERS DETAILS:</b>			
<b>FULL NAME:</b>			
<b>PURCHASING DATE:</b>		<b>PHONE::</b>	
<b>SERIAL NUMBER:</b>			
<b>LOT :</b>			
<b>DEALER DETAILS</b>			
<b>FULL NAME:</b>			
<b>DDRESS:</b>			
<b>PHONE:</b>			
<b>SIGN AND SIGNATURE</b>			

### INCIDENT REPORTING:

Any serious incident that occurs in relation to use of the device must be reported to the manufacturer and a competent authority of the Member State in which the user and / or patient resides.

MOBIAC S.A  
ΚΑΘΙΑΝΑ ΑΚΡΩΤΗΡΙΟΥ, ΧΑΝΙΑ - ΚΡΗΤΗ  
KATHIANA AKROTIRIOU CHANIA - CRETE  
T: +30 28210 63222 | WEB: www.mobiakcare.com